

ENT


Middle Ear Instruments	60
Picks and Scoops	62
Hooks and Probes	63
Retractors	64
Suction Tubes	65
Speculum	66
Aural Forceps	68
Artery Forceps	70
Tonsil, Turbinate and Septum Forceps	71
Dressing Forceps	72
Other Forceps	73
Mouth Gags	74
Tongue Plates and Depressors	75
Punches	76
Rongeurs	80
Scissors	82
Gouges, Knives and Elevators	84
Curettes and Myringotomes	85
Elevators and Dissectors	86
Miscellaneous	87

SURGEON'S FIRST CHOICE

Middle Ear Instruments

Ironside Instruments


160mm long, black finish



Description
Knife 1, vertical incision
Knife 2, horizontal incision
Elevator
Needle point
Hook, superior elevation, angle 45°
Hook, inferior elevation, angle 45°
Hook, anterior elevation, angle 90° right
Hook, anterior elevation, angle 90° left
Measuring device

Rosen Instruments


160mm long, black finish



Description
Curette, tympanic ring, small
Curette, tympanic ring, large
Probe, footplate movement test, fine
Mobiliser, slight curve
Mobiliser, full curve
Scoop, size 000
Scoop, size 00
Scoop, size 0
Scoop, size 1
Fenestrator, oval window

Schuknecht Instruments


160mm long, black finish



Description
Hook upward, 0.4mm
Hook downward, 0.4mm
Hook upward, 0.7mm
Hook downward, 0.7mm
Elevator, oval
Hook, right-angled upward
Pick, straight
Spatula

Beale Instruments


165mm long, satin finish



Description
Tympanoplasty probe, fine blunt end
Raspatory, slight curve
Raspatory, slight curve, small
Rugine, double ended, sharp
Rugine, double ended
Spatula, graft smoothing

Beale Instruments


170mm long, satin finish



Description
Knife, Stapedius tendon cutting
Raspatory, mucous removal
Osteotome, footplate
Perforator, elevation fixed footplate
Hook, removal of footplate, angled downward
Hook, removal of footplate, angled upward
Hook, bone removal from oval window, angled left
Applicator, vein graft for right ear
Applicator, vein graft for left ear
Curette, removal bony rim
Curette

Hough Instruments


165mm long, stainless steel



Description
Excavator, "whirlybird", right, large
Excavator, "whirlybird", left, large
Excavator, "whirlybird", right, small
Excavator, "whirlybird", left, small
Pick, stapedectomy footplate, angled 90°, 0.33mm
Pick, stapedectomy footplate, angled 45°, 0.5mm
Pick, stapedectomy footplate, angled 15°, 0.5mm
Hook, "hoe", excavator, angled 90°
Hook, "hoe", excavator, angled 45°
Pick, curved

Cawthorne Hooks

165mm long, black finish



Description
Slight curve
Medium curve
Full curve
Double curve
90° angle, 0.5mm
90° angle, 1.0mm
45° angle, 0.5mm
Retrograde
Straight

Protect your delicate ENT instruments with our custom baskets

Manufactured at our workshop in Essex



Contact us with your requirements

info@syze.co

Picks and Scoops

Musgrove Cerumen Pick and Wool Carrier

Double ended



Description

178mm long

Lempert Scoop

Oval, 200mm long



Description

2mm x 2.5mm

2mm x 3.5mm

2mm x 4mm

2mm x 4.5mm

3mm x 5.5mm

3mm x 6mm

Formby Cerumen Scoop and Hook

Double ended



Description

178mm long

House Curette

Double ended, 180mm long



Description

Straight

Angled

Hovell Loop

Fine, for removing hard wax



Description

152mm long

Balance Scoop and Curette



Description

230mm long

Hartmann Hook

Sharp, straight shaft



Description

140mm long

Gardiner Brown Hook and Imray Loop

Double ended



Description

127mm long

Cawthorne Cerumen Hook



Description

143mm long

St Bart's Hook

Blunt, straight handle



Description

178mm long

Tracheal Hook

Sharp



Description

165mm long

McIndoe Skin Hook

Fine, with flat handle



Description

191mm long

Tracheal Single Hook Retractor

Sharp



Description

152mm long

Cawthorne Probe

Silver, flattened end



Description

140mm long

Tracheal Double Hook

Blunt



Description

105mm long

Jobson Horne Probe

Double ended, serrated tip and smooth ring



Description


140mm long

178mm long

Retractors

Mollison Retractor

Self retaining, curved




Description

4-4 teeth, adult, 165mm long

2-2 teeth, child, 133mm long

Wallace Modified Plester Retractor




Description

2 x 2 teeth, 140mm long

Plester-Jansen Retractor

Solid blade one side,
140mm long



Description


Left

Right

3 x 3 Teeth

Shah Intra Aural Retractor

Satin finish
152mm long



Description

4-4 teeth, adult, 165mm long

Third blade

Intra aural blade

Post aural blade

Hurd Pillar Retractor and Tonsil Dissector



Description

197mm long

Mollison Anterior Pillar Retractor and Tonsil Enucleator

Semi sharp




Description

210mm long

SURGEON'S FIRST CHOICE

Ipswich Pattern Sucker



Description

267mm long

Lempert Suction and Irrigation Tubes

With tucker valve and stilette for cleaning stainless steel, 203mm




Description

6F.G. 2mm

9F.G. 3mm

Formby Suction Tube

Angled, conical, tip size 10 F.G. stainless steel



Description

210mm long

Magill Suction Tubes

Set of four with conical ends, one curved, three angled, handle, stainless steel, 210mm



Description

Various sizes

3.0mm ø

3.5mm ø

4.0mm ø

Barron Suction Tube

Small size, 4 F.G. Not for use with Zoellner ends



Description

229mm long

Zoellner Suction Tube

With stilette, stainless steel



Description

170mm long, 2mm ø

House Suction and Irrigation Tubes

Stainless steel, 152mm long




Description

8 F.G. x 6 F.G.

10 F.G. x 7 F.G.


12 F.G. x 8 F.G.

Rowbotham Suction Tube



Description

286mm long



Speculum

Rosen Slotted Speculum

Set of four, round ends, black finish
38mm long


Description
Set of 4
4mm
5mm
6mm
7mm



Tumarkin Speculum

Set of eight, tapered slot, black finish
50mm long

Description
Set of 8
4.5mm ø
6.5mm ø
7.0mm ø
7.5mm ø
8.0mm ø
8.5mm ø
9.0mm ø
11.0mm ø



Rosen Aural Speculum

Round, 40mm long

Description
3mm
4mm
5mm
6mm
7mm



Heath Speculum

Set of four, black finish
50mm long

Description
Set of 4
4.0mm ø
5.0mm ø
6.0mm ø
7.0mm ø



Shea Speculum

Set of five, oblique ends, black finish, 38mm long

Description
Set of 5
3mm
4mm
5mm
6mm
7mm




Peters Seigle Speculum

Description
Peters ear speculum black chrome finish set of 5, sizes 0 to 4
2.2mm x 3.2mm
2.8mm x 4.2mm
4.3mm x 5.8mm
5.2mm x 7.6mm
6.0mm x 7.6mm
6mm
7mm
8mm
9mm
10mm



St Clair Thompson Nasal Speculum

Description
38mm long
51mm long
64mm long
89mm long



Killian Nasal Speculum

With Stop

Description
51mm blade
64mm blade
76mm blade
89mm blade




Without Stop

Description
76mm blade

Thudichum Nasal Speculum

Set of seven, sizes 1-7,
70mm long

Description
Standard
Insulated




Standard

Description
Size 1 - 12mm blade
Size 2 - 14mm blade
Size 3 - 16mm blade
Size 4 - 18mm blade
Size 5 - 20mm blade
Size 6 - 22mm blade
Size 7 - 24mm blade
Size 8 - 27mm blade

Insulated

Description
Size 1 - 12mm blade
Size 2 - 14mm blade
Size 3 - 16mm blade
Size 4 - 18mm blade
Size 5 - 20mm blade
Size 6 - 22mm blade
Size 7 - 24mm blade
Size 8 - 27mm blade

Lempert Endaural Speculum



Description
133mm long

Pilcher Nasal Speculum



Description
152mm

Cottle Nasal Speculum

127mm long



Description
20mm blade
32mm blade
41mm blade
60mm blade

**Consistent quality across
all our instruments**



Did you know...

Surgical Holdings have patterns and drawings to ensure their instruments are consistent in design and finish for end users.

Get in touch to discuss your instrument requirements at info@syze.co

Aural Forceps



Beale Aural Forceps

Micro, black finish, triangular flat jaws, serrated



Description
70mm to shoulder

Hough Aural Forceps

Micro, extra fine oval cup jaws, black finish



Description
70mm to shoulder

Cawthorne Aural Forceps

Micro, black finish, extra fine
70mm to shoulder



Description
Smooth jaws
Serrated jaws

Ormerod Aural Forceps

Micro, extra fine, black finish
70mm to shoulder



Description
Serrated, angled upwards
1mm cup jaw
Elongated cup jaw, 0.5 x 4mm

Hartmann Aural Forceps

Micro, black finish
75mm to shoulder



Description
Serrated jaws
Tenaculum jaws
Cup shape jaws
Curved left, serrated
Curved right, serrated
Angled up, serrated
Serrated jaws, satin finish
Pre-sterile - single use only

Shea Aural Forceps


Micro, oval cup jaw, black finish
70mm to shoulder



Description
Standard
Curved to left
Curved to right
Curved upward

Hartmann Aural Forceps

Standard, satin finish
75mm to shoulder




Description

Serrated jaws
Tenaculum jaws
Cup shape jaws
Grunwald punch action

Heath Aural Forceps

Light, with lifting discs, serrated jaws




Description

83mm to shoulder

Henckel Aural Forceps

Fine fenestrated scoop blades




Description

75mm to shoulder

Struycken Cutting Forceps

75mm to shoulder



Description

Fenestrated under jaw

SURGEON'S FIRST CHOICE

Did you know...

We can repair delicate ENT instruments

Our technicians manufacture instruments from scratch, so they have the skills and experience needed to repair even the most delicate instruments.


Get in touch to book a repair

info@syze.co

Artery Forceps

Birkett Tonsil Artery Forceps

Fine, 178mm long




Description

- Straight
- Curved

Negus Tonsil Artery Forceps


180mm long



Description

- Large curve
- Small curve

Denis Browne Tonsil Artery Forceps




Description

216mm long

Sawtell Tonsil Artery Forceps

178mm long




Description

- Large curve
- Small curve

Fickling Artery Forceps

190mm long




Description

- 1-2 teeth
- Serrated

Scott Tonsil Artery Forceps

Curved to side, screw joint




Description

180mm long

Tonsil, Turbinate and Septum Forceps

Denis Browne Tonsil Holding Forceps

Angled shanks, no rack




Description

200mm long

Jansen-Middleton Nasal Septum Forceps

Compound action
178mm long




Description

5mm jaw

Irwin Moore Nasal Turbinate Forceps


191mm long



Description

5mm jaw

Walsham Septum Forceps




Description

Straight, 215mm long

Right, 215mm long

Left, 215mm long

Luc Nasal Turbinate Forceps



Heart shaped, 178mm long

Description

Small 6mm ø

Medium 7.5mm ø

Large 9.0mm ø

Oval, 191mm long

Description


Small 6mm ø

Medium 7.5mm ø

Large 9.0mm ø

Hartmann Conchotomes

Circular cutting
102mm to shoulder



Description

5mm ø

7mm ø


9mm ø

11mm ø

Dressing Forceps

Hartmann Dressing Forceps

Serrated jaw, screw joint



Description
108mm to shoulder

Wilde Dressing Forceps


Angular, serrated jaws



Description
102mm to shoulder

Tilley Dressing Forceps


Straight, fine, serrated ends



Description
2mm wide, 76mm to shoulder

Irwin Moore Thyroid Clamp

Thin jaws, angled to side



Description
165mm long

Secto® Dissectors

Designed to provide a gentle, atraumatic way to separate tissue during surgical procedures. These sterile, cotton-constructed dissectors are x-ray detectable and available in four patterns: Cherry, Triangle, Peanut and Kittner.



Secto® Cherry Dissectors



Secto® Triangle Dissectors



Secto® Peanut Dissectors



Secto® Kittner Dissectors

Sterile, cotton, non-linting, x-ray detectable gauzes with identification tape and string for use during tonsillectomies.




Secto® Tonsil Sponges



Contact us with your requirements info@syze.co

Trousseau Tracheal Dilator



Description

Adult 140mm long

Child 125mm long

Wilson Bone Nibbling Forceps

Fine, curved, satin finish
152mm long




Description

Downward cutting - 3mm jaws

Upward cutting - 3mm jaws

Cawthorne Nibbling Forceps


Narrow, 3mm jaws, wide handle



Description

159mm long

Denis Browne Collar Crushing Forceps



Description

140mm long

Lempert Nibbling Forceps

152mm long



Description


Straight, 2.5mm jaws

Angled, 3.0mm jaws

SURGEON'S FIRST CHOICE


Mouth Gags

Doyen Mouth Gag



Description
Adult
Child


Fergusson Ackland Mouth Gag



Description
Adult, 210mm
Child, 155mm


Boyle-Davis Mouth Gag

Frame only



Description
Adult
Child

Mason Mouth Gag



Description
Adult, 210mm

Draffin Bipod



Adult


Description
Adult , 4 hole x 480mm long half circle end for use with a drape
Adult , 4 hole x 480mm long with round ends to be used with Magauran plate
Adult , 4 hole x 480mm long with round ends to be used with Magauran plate

Child

Description
Child, 3 hole x 425mm long half circle end for use with a drape
Child, 3 hole x 425mm long with round ends to be used with Magauran plate
Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern
Child, 3 hole x 425mm long with round ends to be used with Magauran plate
Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern
Magauran plate 480mm - for use with Draffin bipods

Kilner Dott Mouth Gag


Frame only



Description
Adult
Set, including frame, small, medium, and large blades


Magauran Plate

Description
480mm - for use with Draffin bipods



Tongue Plates and Depressors

Doughty Tongue Plate



Plain

Description
65mm - 1.6mm-1.9mm slot
75mm - 1.6mm-1.9mm slot
90mm - 4mm-4.5mm slot
100mm - 4mm-4.5mm slot
115mm - 4mm-4.5mm slot

Slotted

Description
65mm - 1.6mm-1.9mm slot
75mm - 1.6mm-1.9mm slot
90mm - 4mm-4.5mm slot
100mm - 4mm-4.5mm slot
115mm - 4mm-4.5mm slot

With anaesthetic tube

Description
65mm - 1.6mm-1.9mm slot
75mm - 1.6mm-1.9mm slot
90mm - 4mm-4.5mm slot
100mm - 4mm-4.5mm slot
115mm - 4mm-4.5mm slot

Kilner Dott Tongue Plate

With endotracheal tube groove



Description
Small
Medium
Large
Extra small


Tongue Depressor

Flat stainless steel



Description
184mm long

Lack Tongue Depressor



Description
Child, 13mm wide
Adolescent, 19mm wide
Adult, 22mm wide

Punches


Nasal Cutting Forceps



Myles


Through cutting, 120mm to shoulder

Description
12.5 x 3.5mm jaw



Seymour Jones Antrum Punch Forceps

159mm long




Description
4mm jaw

Grunwald-Henke

Through cutting, 120mm to shoulder


Description
7 x 2mm jaw
9 x 2.5mm jaw



Struycken


Through cutting, 120mm to shoulder

Description
11 x 1mm jaw



Backbiter


180mm long



Description
7 x 3mm jaw

Grunwald Nasal Punch Forceps


108mm to shoulder



Description
A fenestrated jaw, 3.5mm
Fenestrated lower jaw, narrow, 2.5mm
Fenestrated and tapered lower jaw, 2-3mm

Backbiter


200mm long



Description
7 x 2mm jaw

Stamberger Rhinoforcep Antrium Punch

7mm x 3mm, 130mm long




Description

Left back cut
Right back cut

Tilley Henckel Ethmoid Punch

Graduated shank,
114mm to shoulder




Description

4.5mm
5.5mm
6.5mm
8mm

Ostrom Antrum Punch

Forcep, 178mm to shoulder




Description

3mm jaw

Mackay Modified Stamberger Punch

200mm long




Description

7 x 2mm jaw

SURGEON'S FIRST CHOICE

Paterson Biopsy Forceps

Elongated cup, 3mm jaw



Description


305mm long
355mm long
406mm long
508mm long

Punches

Chevalier Jackson Endoscope Forcep

With crocodile action and serrated jaw

Description
Effective shaft length 300mm
Effective shaft length 400mm
Effective shaft length 500mm
Effective shaft length 600mm



Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and curved cup

Description
Effective shaft length 300mm
Effective shaft length 400mm
Effective shaft length 500mm
Effective shaft length 600mm



Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and straight cup


Description
Effective shaft length 300mm
Effective shaft length 400mm
Effective shaft length 500mm
Effective shaft length 600mm



St Bartholomew's Ethmoid Punch


Scoop blades, 3mm bite, 114mm to shoulder

Description
Straight
Upward cutting
Downward cutting




Ferris Smith Kerrison Punch

Heavy, 5mm bite,
229mm long



Citelli Antrum Punch

76mm to shoulder



Description


Upward cutting
Downward cutting

Description

3mm jaw
5mm jaw

Hajek Antrum Punch

197mm long



Description

5mm jaw

Hajek Sphenoidal Punch

146mm to shoulder



Upward cutting

Description

4mm bite
5mm bite

Downward cutting

Description

4mm bite
5mm bite

Rongeurs

Beck Rongeur

Straight, for intervertebral disc



Description
2mm bite, 130mm shaft

Pennybacker Rongeur


120mm shaft



Description
3mm bite
4mm bite

Cushing Rongeur


130mm shaft



Description
2mm bite
3mm bite
4mm bite
5mm bite
6mm bite

St Bart's Rongeur

114mm shaft



Description
3mm bite

Cushing Rongeur

130mm shaft



150° upward cutting
2mm bite
3mm bite
4mm bite
5mm bite


150° downward cutting
2mm bite
3mm bite
4mm bite
5mm bite

St Bart's Rongeur

140°, 3mm, 114mm shaft



Description
Upward cutting
Downward cutting



Mackay Rongeur

140°, straight, 114mm shaft




Description

5.5mm

Weil-Blakesley Rongeur

110mm shaft



140° upward cutting

Description

3.6mm bite
4.2mm bite
4.8mm bite


140° downward cutting

Description

3.6mm bite
4.2mm bite
4.8mm bite

Mackay Rongeur

With suction tube



Description

3.5mm bite, 114mm shaft

Weil-Blakesley Rongeur

120mm shaft




Description

3.0mm bite
3.6mm bite
4.2mm bite
4.8mm bite
5.8mm bite
6.6mm bite

Mackay Rongeur

With suction tube, 140° upward cutting



Description

3.5mm bite, 114mm shaft

Scissors

Beckmann Turbinectomy Scissors

Curved




Description

165mm long

Fomon Scissors


Straight, serrated blades, angled handles



Description

145mm long

Scott Stevenson Tonsil Scissors




Description

216mm long

Cottle Scissors

Straight blades, angled handles, blunt/blunt



Description

165mm long

Wood Tonsil Scissors

Sharpened outer edges




Description

178mm long

Cottle Bulldog Scissors

Blunt/blunt, 115mm long




Description

115mm long

Wheeler Tonsil Scissors

Triangular form blades, serrated shanks




Description

184mm long

Joseph Scissors

Sharp/sharp, 150mm long




Description

Straight

Curved

Tilley Tonsil Scissors

Curved




Description

190mm long

Heymann Turbinectomy Scissors


Angular



Description

197mm long

Wilson Tonsil Scissors



Description

203mm long

Cawthorne Scissors

Micro, extra fine, 70mm to shoulder



Description


Left

Right

Straight



Thyroid Scissors



Description

152mm long

Ormorod Aural Scissors

Micro, extra fine, curved upwards



Description

70mm to shoulder



Knapp Scissors

Slight curve, blunt points, 114mm long



Description

Satin finish

Black finish

Portmann Aural Scissors

Micro, satin finish 70mm to shoulder



Description

Curved to right

Curved to left




SURGEON'S FIRST CHOICE

Gouges, Knives and Elevators

Ballenger Nasal Cartilage Knife

Revolving blade



Straight, 203mm long

Description
6mm wide
8mm wide

Cranked, 197mm long

Description
6mm wide

Guy's Hospital Nasal Gouge

Shallow groove and flattened back, 127mm long



Description
4mm

Tilley Nasal Gouge

V shaped, 127mm long



Description
4mm

Sickle Knife



Description
188mm long
160mm long

Cottle Elevator

Double ended, black chrome finish
200mm long



Jenkins Gouge

140mm long



Description
1mm
2mm
3mm
4mm
5mm
6mm
7mm

Description
8mm
9mm
10mm
11mm
12mm
13mm
14mm

Hill Elevator

Double ended, curved to side, rugine end



Description
5mm wide, 222mm long

Killian Nasal Gouge

127mm long



Description
4mm

Tumarkin Atticotomy Elevator

Double ended




Description
191mm long

Curettes and Myringotomes

Buck Ear Curette

Angled, blunt, 146mm long



Description
2mm
2.5mm
3mm
3.5mm
4mm

Nasal Curette


Straight, 185mm long



Description
5 x 2mm

Beckmann Adenoid Curette


216mm long



Description
8mm blade
10mm blade
12mm blade
14mm blade
16mm blade
18mm blade

Agnew (Poltzer) Myringotome

190mm long



Description
Shaft length 100mm


Trautmann Myringotome

178mm long



St Clair Thomson Adenoid Curette

With cage, 235mm long




Description
8mm
10mm
12mm
14mm
16mm
18mm

Description
Upward cutting
Downward cutting

Elevators and Dissectors

Freer Instruments

Single ended, 190mm long



Knife

Description
Curved to side
Rounded edge, 8mm blade
Flat rounded blade
Rounded blade, curved on flat


Elevator

Description
Dull edged
Sharp edged

Chisel

Freer Separator and Elevator

Double ended



Description
Blunt/sharp, 180mm long
Sharp, knife curved to side, 180mm long
Blunt/sharp, 200mm long
Blunt/sharp, 188mm long

St Clair Thompson Dissector

159mm



Description
Sharp
Blunt

Negus Aspiring Dissector

With 6mm olive, 210mm long



Description
Insulated

Syme Tonsil Dissector

Double ended, one end blunt, one semi sharp



Description
170mm long

Gwynne-Evans Tonsil Dissector

Double ended, with serrated end and blunt end 197mm long



Description
6mm wide
8mm wide

Woakes Eustachian Catheters

Silver, Luer, 140mm




Description
7 F.G.
8 F.G.
9 F.G.

Description
10 F.G.
11 F.G.

Laryngeal Mirror

Mirror only, 130mm long



Description
6mm ø
8mm ø
10mm ø
12mm ø
14mm ø
16mm ø
18mm ø

Description
20mm ø
22mm ø
24mm ø
26mm ø
28mm ø
30mm ø

Kramer Eustachian Catheters

Silver, Luer, 140mm




Description
7 F.G.
8 F.G.
9 F.G.

Description
10 F.G.
11 F.G.

Laryngeal Mirror

Complete with metal handle, 260mm long



Myles Naso-Antral Perforator

171mm long



Description
16 F.G.
19 F.G.

Description
6mm ø
8mm ø
10mm ø
12mm ø
14mm ø
16mm ø
18mm ø

Description
20mm ø
22mm ø
24mm ø
26mm ø
28mm ø
30mm ø


Mitchell Trimmer



Description
165mm

Tilley Antrum Burrs

171mm long



Description
6mm ø
7mm ø
8mm ø

Miscellaneous

Lake Seeker

Double ended, for vestibule



Description

184mm long

Negus Pusher


"Knot tier", metal handle



Description

203mm long


Lempert Rugine



Description

Narrow, 156mm long
Broad, 171mm long

Stanford Cade Insufflator



Description

175mm long

Irwin-Moore Tonsil Needle

203mm long



Description

Curved to right
Curved to left

Tubing Mount

For use with catheters, Luer



Description

25mm