

Gynaecology

Clamps	90
Forceps	90
Biopsy Forceps	94
Cervical Forceps	97
Obstetric Forceps	98
Speculum	99
Retractors	101
Curettes	103
Catheters	104
Scissors	106
Miscellaneous	107
Dilators and Sounds	108




SYRZA

SURGEON'S FIRST CHOICE

Clamps and Forceps

Parametrium Clamp Forceps

Hysterectomy clamp, atraumatic jaws, two gold bows



Straight

Description
203mm long
233mm long
254mm long
304mm long
360mm long

Slight curve

Description
203mm long
233mm long
254mm long
304mm long
360mm long

Strong curve

Description
203mm long
233mm long
254mm long
360mm long

Angled

Description
203mm long
233mm long
254mm long
304mm long
360mm long

Spencer-Robarts Clamp


For endometrial and pedunculated fibroid polyps, with ratchet jaw-locking mechanism, and 50° jaw opening



Description
Shaft 235mm long

Museux Vulsellum Forceps

Matt finish




Description
240mm long



Schroeder Uterine Vulsellum Forcep


1-1 teeth



Description
229mm long

Teale Uterine Vulsellum Forcep

3-4 teeth, 230mm long



Description
Curved
Straight

Bonney Hysterectomy Forceps

Angled on flat, 1 - 2 teeth




Description

197mm long

Maingot Hysterectomy Forceps

203mm long




Description

Curved on flat, 1 - 2 teeth

Straight

Berkeley-Bonney Hysterectomy Forceps


Curved on flat, 1 - 2 teeth



Description

203mm long

Faures Peritoneum Forceps



205mm long

Description

Slight curve

203mm long,
longitudinal serrations


Description

Curved, toothed

Curved

Gwilliam Hysterectomy Forceps

Longitudinal serrations
203mm long




Description

Straight,

Curved on flat, 1-2 teeth

Heaney Hysterectomy Forceps

Curved on flat, 203mm long



Description

One tooth


Double tooth

SURGEON'S FIRST CHOICE

Forceps

Bonney Uterine Polypus Forcep

Box joint




Description

240mm long

Somer Uterus Holding Forcep

Matt finish




Description

230mm long

254mm long

Noto Uterine Polypus Forceps

Matt finish



Description

240mm long

McClintock Ovum Forcep

Curved, screw joint




Description

241mm long

Vant Uterine Polypus Forcep

Serrated




Description

240mm long

Greenhalgh Ovum Forcep

Screw joint




Description

240mm long

Dartigues Uterine Forceps

Matt finish




Description

270mm long

Keily Ovum Forceps

Matt finish



Description

320mm long

Doyen Ovarian Forceps

Matt finish



Description

210mm long

Lowa Tenaculum Forceps

Matt finish, 5-6 teeth




Description

260mm long



Wertheim Parametrium Forcep

Curved on flat, screw joint




Description

254mm long

Skene Tenaculum Forceps

Bright finish




Description

240mm long

Green-Armytage Haemostasis Forcep

Screw joint, 15mm jaw, 210mm long



Description

Straight
Angled

Pratt Tenaculum Forceps

Matt finish



Description

280mm long

Kogan Endocervical Forceps

Matt finish, 275mm long




Description

3mm jaw
5mm jaw
7mm jaw

Coil Removing Forcep

With serrated rounded jaws
4mm wide




Description

160mm to shoulder

Forceps and Biopsy Forceps

Willet Martel Forceps

Matt finish




Description

190mm long



Littlewood Tissue Forceps

2-3 teeth, matt finish




Description

190mm long

Eppendorf Biopsy Punch Forcep

200mm to shoulder




Bite shown actual size

Description

3mm x 8mm bite

Kevorkian Biopsy Punch Forcep

255mm to shoulder




Bite shown actual size

Description

3mm x 8mm bite

Tischler Morgan Biopsy Forcep

210mm to shoulder




Bite shown actual size

Description

3mm x 7mm bite

Tischler Morgan Biopsy Punch Forcep

295mm to shoulder



Bite shown actual size


Description


2mm x 7mm bite

Biopsy Forceps

Baby Tischler Biopsy Forcep

210mm to shoulder




 Bite shown actual size


Description

2mm x 5mm bite

Baby Tischler Biopsy Punch Forcep

295mm to shoulder




 Bite shown actual size


Description

2mm x 5mm bite

Van Doren Biopsy Forceps


Bright finish, 260mm overall length



 Bite shown actual size


Description

5mm x 9mm bite



Van Doren Cervical Biopsy Punch Forcep

With curved basket jaw, angled shanks




Description

240mm long

Wittner Cervical Biopsy Punch Forcep

240mm long shaft



Description

Jaws curved to side


Jaws straight to side

SURGEON'S FIRST CHOICE

Obstetric Forceps

Barnes Obstetric Forcep

With Simpson handles and Neville axis traction rod




Description

381mm long

Haig Ferguson Obstetric Forcep

Tape slots in blades and axis traction rod




Description

356mm long

Kieland Obstetric Forcep

Matt finish




Description

410mm long

Simpson Braun Obstetric Forcep

Matt finish




Description

380mm long

Wrigley Obstetric Forcep

Short model for caesarian section




Description

279mm long

Auvard Vaginal Speculum


Weighted



Description

Small blade 38 x 90mm
Large blade 51 x 90mm

Collin Vaginal Speculum



Description

Small - blades 96mm x 25mm
Medium - blades 108mm x 30mm
Large - blades 108mm x 35mm


Graves Vaginal Speculum

Bright finish



Description

115 x 35mm



Nott Speculum


Stainless steel

Description

Blade 110mm long, 30mm wide

Cusco Vaginal Speculum

Stainless steel



Description

Extra small - 98mm x 16mm
Small - 105mm x 25mm
Medium - 107mm x 29mm
Large - 115mm x 35mm

Insulated

Description

Small - 105mm x 25mm
Medium - 107mm x 29mm
Large - 115mm x 35mm
Extra large - 135mm x 38mm

With Suction Tube

Description

Small - 105mm x 25mm
Medium - 107mm x 29mm
Large - 115mm x 35mm
Extra large - 135mm x 38mm

Black Finish


Description

Small - 105mm x 25mm
Medium - 107mm x 29mm
Large - 115mm x 35mm
Extra large - 135mm x 38mm

Speculum

Pederson Vaginal Speculum

Stainless steel




Description

120mm x 25mm wide

90mm x 20mm wide

Winterton Vaginal Speculum

Stainless steel




Description

Blade 115mm long, 25mm wide

Modified Speculum

Stainless steel, for elderly patients




Description

Blade 180mm long, 18mm wide

Sims Duckbill Speculum

Double ended



Description


Small - 27mm x 29mm x 130mm long

Medium - 30mm x 33mm x 158mm long

Large - 33mm x 38mm x 184mm long

Sims Duckbill Speculum

Double ended,
insulated



Description

Small - 27mm x 29mm x 130mm long

Medium - 30mm x 33mm x 158mm long

Large - 33mm x 38mm x 184mm long



Bovie Smoke Shark II

Optimal Control, Better Protection

Smoke Shark is essential in procedures where HPV may be present and could be inhaled by either patient or clinician


- Protects OR staff from the risk of exposure to Biological, Physical and Chemical Hazards
- Lightweight, compact and ultra-quiet
- Easy to set up and operate
- Extensive life filter (SF35) offers up to 35 hours performance with 3 flow settings



info@syze.co

Breisky Vaginal Retractor


Bright finish, matt handle
320mm long



Description

25mm blade


Doyen Retractor



Description

Blade 100mm x 30mm


Doyen Rectangular Retractor



Description

Small - 48 x 60mm
Medium - 48 x 89mm
Large - 48 x 124mm

St Mark's Hospital Pattern Retractor



With lip

Description


Blade 125mm x 60mm
Blade 175mm x 60mm

Without lip

Description

Blade 125mm x 60mm
Blade 175mm x 60mm


Doyen Retractor



Description

Blade 95mm x 27mm

Wertheim Retractor




Description

Blade 95mm x 38mm

Retractors

Landon Vaginal Retractor



Description

Blade 89mm x 25mm

Lateral Vaginal Retractor

Blades bright finish, matt handle



Description

60 x 20mm blade

75 x 22mm blade

110 x 23mm blade

AVAILABLE INSULATED
ADD SUFFIX 'INS'

Whitlow Retractor

Vaginal retractor, for use with supplied blades. Points of retraction move around almost 360°. Supplied complete with sterilising basket and three sets of blades.

Frame



225mm at widest point

Blades

1 x anterior and 1 x posterior 90°, 25mm x 100mm

2 x lateral 90°, 25mm x 70mm

2 x lateral 75°, 25mm x 70mm



Kevorkian Biopsy Curette

313mm long

Description
With basket
Without basket



Novak Curette

Description
230mm long



Recamier Uterine Curette

Matt finish, 310mm long

Description
Blunt
Sharp



Rheinstadter Flushing Curette

290mm long

Sharp

Description
Sharp, 8mm
Sharp, 14mm

Blunt

Description
Blunt, 8mm
Blunt, 14mm



SURGEON'S FIRST CHOICE

Sims Uterine Curette

Single ended, 255mm long, matt finish

Sharp

Description
Small 6mm
Medium 10mm
Large 13mm

Blunt

Description
Small 6mm
Medium 10mm
Large 13mm



Sims Uterine Curette

Double ended, 250mm long

Sharp

Description
Small 7 x 5mm
Medium 10 x 7mm
Large 13 x 10mm

Blunt

Description
Small 7 x 5mm
Medium 10 x 7mm
Large 13 x 10mm



Catheters

Rotunda Hospital Pattern Catheter

With solid end, chromium plated, 203mm long



Description

8 E.G.

10 E.G.

Catheter

Oval, solid end, chromium plated



Description

165mm long

Catheter Introducer

457mm long



Description

Straight, 2mm ϕ

Curved, 2mm ϕ

Catheter Introducer



Description

Fine, 395mm long

Catheter Introducer

Bulbous end, with eye in handle



Description

Curved, 2mm x 410mm long

Catheter Introducer

With eye in handle, curved



Description

2mm x 420mm long

Catheter Introducer

Fixed handle



Description

445mm long

Catheter Introducer



Description

400mm long

Catheter Introducer

410mm long



Description

410mm long

Scissors

Barnes Episiotomy Scissors




Description

127mm long

Braun Stadler Episiotomy Scissors

Matt finish



Description

140mm long

222mm long

Vant Episiotomy Scissors



Description

76mm blades, 241mm long

Busch Umbilical Scissors

Matt finish



Description

160mm long

Umbilical Scissors

American pattern




Description

102mm long

Simms Uterine Scissors

Curved, tungsten carbide, with two gold bows




Description

200mm long



Dubois Uterine Scissors

Matt finish



Description

270mm long

Uterine Scissors

Curved on flat



Description

203mm long

Uterine Scissors

Angled on flat





Description

230mm long

Parametrium Scissors

225mm long



Tungsten carbide, standard, with two gold bows

Description
Slight curve
Strong curve
Acute curve


Tungsten carbide, supercut, with one gold and one black bow

Description
Slight curve
Strong curve
Acute curve

Miscellaneous

Bonney Myomectomy Clamp

With 51mm jaws, screw joint




Description

267mm long

Sims Vaginal Depressor

Double ended



Description

260mm long

SURGEON'S FIRST CHOICE

Spackman Uterine Cannula

Chrome plated, 360mm long



Description

- 12 F.G. with rubber cone
- 14 F.G. with metal cone
- 14 F.G. with rubber cone
- 18 F.G. with metal cone
- 18 F.G. with rubber cone

Doyen Myoma Screw Large



Description

178mm long

Doyen Myoma Screw Small



Description

178mm long

Dilators and Sounds

Galabin Uterine Sound

Calibrated 4-24, bright finish



Description

305mm long

Sims Uterine Sound

Calibrated 4-24, bright finish



Description

320mm long

Umbilical Dilator



Description

95mm long

Hegar Cervix Dilator

Single ended, 185mm long



Description

1mm ø
2mm ø
3mm ø
4mm ø
5mm ø
6mm ø
7mm ø
8mm ø
9mm ø

Description

10mm ø
11mm ø
12mm ø
13mm ø
14mm ø
15mm ø
16mm ø
17mm ø

Hawkin Ambler Cervix Dilators

Conical, 171mm long



Description


Set of 12
Size 3 - 6mm
Size 4 - 7mm
Size 5 - 8mm
Size 6 - 9mm
Size 7 - 10mm
Size 8 - 11mm

Description

Size 9 - 12mm
Size 10 - 13mm
Size 11 - 14mm
Size 12 - 15mm
Size 13 - 16mm
Size 14 - 17mm

Hegar Cervix Dilator

Double ended, 200mm long



Description

Set of 8: from sizes 1 and 2mm to 15 and 16mm
1 and 2mm ø
3 and 4mm ø
5 and 6mm ø
7 and 8mm ø

Description

9 and 10mm ø
11 and 12mm ø
13 and 14mm ø
15 and 16mm ø
17 and 18mm ø