



Neuro/ Spinal Surgery

Forceps	115
Scissors	117
Dissectors and Probes	118
Hooks and Curettes	119
Retractors	120
Elevators	123
Suction Tubes	124
Rongeurs	125

SURGEON'S FIRST CHOICE

Adson Dissecting Forceps

Bayonet, 1-2 teeth, jaw 2.5mm at tip



Description

191mm long

Jefferson Dressing Forceps



Description

Serrated jaws, 178mm long

Adson Dissecting Forceps

Bayonet, fine, serrated jaws, jaw 2.5mm at tip



Description

178mm long

Jefferson Forceps



Description

Fine, 1-2 teeth, 178mm long

Northfield Dressing Forceps

Extra light, serrated jaws



Description


1mm wide, 203mm long

Pennybacker Dissecting Forceps



Description

3-4 teeth, 165mm long



Forceps

Silcock Iris Forceps

102mm long



Description
1-2 teeth
Serrated jaws

Dott Dissecting Forceps

Fine, 1-2 teeth



Description
178mm long

Cairns Dissecting Forceps


Delicate, 2-3 teeth



Description
165mm long
216mm long

Dott Dressing Forceps

Angled



Description
114mm long

Cairns (Dott) Dressing Forceps


Serrated jaws



Description
Medium, 2mm wide, 216mm long
Fine, 1mm wide, 178mm long

Tilley Dressing Forceps


Angled, serrated ends



Description
89mm long

Cairns Scissors

Angled on flat



Description

178mm long

Olivecrona Scissors

Angled, with all micro serrated blade




Description

22mm blade, 150mm shaft

McKissock Scissors

Straight




Description

146mm long

Schmieden (Taylor) Scissors

Angled, probe-pointed under blade




Description

146mm long

Sweet Scissors

Delicate



Description

7mm blade, 144mm shaft

13mm blade, 117mm shaft

SURGEON'S FIRST CHOICE

Dissectors and Probes

Cairns Guide and Director

Malleable



Description

127mm long

Pennybacker Probe Dissector

Double ended



Description

229mm long

Dissector

200mm long



Description

Soft
Soft, with hole

Swedish Dissector

Double ended



Description

Light model, 178mm long
Heavy model, 200mm long

Penfield Dissector



No 1

Broad curve, with 7.5mm sharp round spoon, 184mm long

Watson Cheyne Dissector

Fine, double-ended, with probe end



No 2

7.5mm sharp round spoon with wax packer, 178mm long



Description

127mm long
191mm long
230mm long
255mm long

No 3

Full curve and wax packer, 197mm long



No 4

Slight curve, 4mm wide, 216mm long



Description

No 1
No 2
No 3
No 4

Hooks

Adson Hook

178mm long



Description
Angled
Straight
Straight, blunt

Curettes

Falconer Scoop Curette


Double ended, angled ends



Description
191mm long

Adson Oliver Neuro Hook

191mm long



Description
Angled downward, blunt
Sharp
Angled upward, blunt

American Spinal Fusion Curettes

Blue acetal handles, 254mm long



Straight

Description	
Size 000	2mm x 4mm
Size 00	2.5mm x 4.5mm
Size 0	2.75mm x 5mm
Size 1	3mm x 5.5mm
Size 2	3.25mm x 6mm
Size 3	4mm x 7mm
Size 4	6mm x 9.5mm
Size 5	7mm x 14mm
Size 6	9mm x 17mm

Angled

Description	
Size 000	2mm x 4mm
Size 00	2.5mm x 4.5mm
Size 0	2.75mm x 5mm
Size 1	3mm x 5.5mm
Size 2	3.25mm x 6mm
Size 3	4mm x 7mm
Size 4	6mm x 9.5mm
Size 5	7mm x 14mm
Size 6	9mm x 17mm

ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER

Cairns Tissue Hook

127mm long



Description
Half round, sharp
Angled, sharp

Oliver Tissue Hook

Blunt




Description
127mm long

SURGEON'S FIRST CHOICE

Retractors

Mastoid Retractor

3 x 3 blunt hooks




Description

114mm long

Henderson Retractor

Self retaining, straight, 4 x 4 blunt pointed blades




Description

178mm long

Jansen (Shea) Retractor

Self retaining, blunt blades




Description

114mm long

Cairns Scalp Retractor

Curved, self retaining, swivel blades




Description

140mm long

Cone Retractor

Self retaining, hinged arms, for use in ventriculography




Description

165mm long

Schnitker Scalp Retractor

Curved, self retaining, blunt swivel blades




Description

140mm long

Jefferson Retractor

Self retaining, curved, blunt pointed blades




Description

140mm long

Cone Laminectomy Retractor

Hinged arms, with cam rack



Description

254mm long

Cone Laminectomy Retractor

Hinged arms


Description

Left, 165mm long

Right, 165mm long

Left, 254mm long

Right, 254mm long




Harvey-Jackson Laminectomy Retractor

Hinged arms, blades 38mm wide x 44mm deep

Description

305mm long




Mayo-Adson (Beckmann) Laminectomy Retractor

Hinged arms

Description

292mm long




Northfield (Cairns) Retractor

Hinged arms, one pair of blades

4 prong : 25mm wide x 16mm deep
and one pair of blades

5 prong: 32mm wide x 25mm deep



Description

292mm long

Meyerding Hemi-Laminectomy Retractor

Description

18cm long 7" blade size 50 x 10mm

18cm long 7" blade size 50 x 20mm

18cm long 7" blade size 60 x 10mm

18cm long 7" blade size 60 x 20mm

18cm long 7" blade size 70 x 10mm


18cm long 7" blade size 70 x 20mm

18cm long blade left 90 x 10mm

18cm long blades left 90 x 20mm

18cm long blades right 90 x 10mm


18cm long blades right 90 x 20mm



Retractors


Cairns-Northfield Laminectomy Retractor

203mm arm length with three pairs of interchangeable swivelling blades, 76mm, 64mm, and 51mm



Hibbs Retractor


203mm long



Description
Blades 10mm x 38mm
Blades 16mm x 54mm

Krayenbuhl Retractor


6mm wide



Description
210mm long
240mm long

Cairns 'Rake' Scalp Retractor

Sharp, 210mm long



Description
3 prong, 16mm wide
4 prong, 19mm wide

Love Retractor


6mm wide



Description
210mm long

Cairns (Cushing) Scalp Retractor

Solid blade



Description
16mm wide, 210mm long

Love Retractor


203mm long



Description
Straight, blade 5mm x 5mm
Shaft angled 90° at handle
Shaft angled 45° at handle

Cairns Strip Retractor

Malleable, 203mm long
8 sizes: 6, 8, 10, 11, 13, 16, 19, 25mm



Description
20 SWG
18 SWG

O'Connell Retractor


Double ended



Description
216mm long

Cobb Spinal Elevator

13mm wide, knurled handle



Cobb Spinal Elevator

With blue acetal handles, 279mm long

Description

10mm wide

13mm wide


19mm wide



ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER

Adson Periosteal Elevator

Serrated grip, 165mm long



straight

Description

7mm wide

curved


Description

5mm wide

7mm wide

McKissock Spade Periosteal Elevator

Straight



Description

197mm long

Cushing Periosteal Elevator

Slightly curved end




Description

197mm long

Pennybacker Elevator

Curved, sharp, 171mm long



Description

6mm wide

SURGEON'S FIRST CHOICE

Suction Tubes

Adson Frazier Suction Tube

Angled, fully insulated, 180mm long




Description

8 F.G.


Dandy Rongeur

Straight



Retractor Suction Tube

With stilette, 220mm long



Description


3mm bite, 178mm shaft

Description

Replacement stilette


Falconer Rongeur

Oblique jaws



Frazier Suction Tubes

Malleable, silver, angled, 180mm long



Description

7 F.G.

8 F.G.

9 F.G.

Description

Upward cutting, 3mm bite, 110mm shaft


Upward cutting, 5mm bite, 110mm shaft

Downward cutting, 3mm bite, 85mm shaft

Downward cutting, 5mm bite, 110mm shaft

Beck Rongeur

Straight, for intervertebral disc




Description

2mm bite, 130mm shaft

Love-Grunwald Rongeur

Cup shaped, jaws 3mm x 10mm, 178mm shaft




Description

- Straight
- Curved upward
- Curved downward

Cushing Rongeur

130mm shaft




Description

- 2mm bite
- 3mm bite
- 4mm bite
- 5mm bite
- 6mm bite

Mackay Rongeur

140°, straight, 114mm shaft




Description

5.5mm

SURGEON'S FIRST CHOICE

Cushing Rongeur

130mm shaft



150° upward cutting

Description

- 2mm bite
- 3mm bite
- 4mm bite
- 5mm bite


150° downward cutting

Description

- 2mm bite
- 3mm bite
- 4mm bite
- 5mm bite

Northfield Rongeur

Curved, upward cutting




Description

4mm bite, 133mm shaft

Rongeurs

Pennybacker Rongeur

120mm shaft




Description

3mm bite
4mm bite

Daniel Rongeur

Miniature, for neuropathology




Description

1.5mm bite, 127mm long

St Bart's Rongeur

114mm shaft




Description

3mm bite

Cairns Rongeur

Fine jaws, angled on flat, serrated jaws, curved handles




Description

4mm bite 152mm long

St Bart's Rongeur

140°, downward cutting, 114mm shaft




Description

3mm

Danish Pattern Rongeur

Angled on flat, fine jaws, curved handles, screw joint




Description

4mm bite, 203mm long

St Bart's Rongeur

140°, upward cutting, 114mm shaft




Description

3mm

Falconer Rongeur

Angled to side, light, screw joint




Description

6mm bite, 203mm long

Wilms Rongeur

Angled on flat, narrow jaws




Description

4mm bite, 203mm long

Pennybacker Rongeur

Straight, fine jaws, 203mm long




Description

4mm bite

Glasgow Pattern Rongeur

Slightly curved on flat, screw joint




Description

4mm bite, 146mm long

Northfield Rongeur

Straight, 10mm bite




Description

One shallow blade, 200mm long

Heanley Rongeur

Long fine straight jaws, screw joint




Description

4mm bite, 178mm long

Stookey Rongeur

Angled to side, fine jaws




Description

3mm bite, 178mm long

Jansen Rongeur

Short curve, screw joint




Description

8mm bite, 152mm long

Trotter Rongeur

Angled to side, strong, screw joint





Description

10mm bite, 203mm long

Rongeurs

Luer-Jansen Rongeur

Light, curved on flat, compound action




Description

4mm bite, 197mm long

Olivecrona Rongeur

Angled to side, compound action, 229mm long





Description

5mm bite

6mm bite

Ruskin Rongeur

Curved on flat, compound action




Description


6mm bite, 184mm long

Leksell Rongeur

Curved to side, compound action




Description
Small, 6mm bite, 229mm long
Large, 8mm bite, 254mm long




Luer Stille Bone Rongeur

Curved, 250mm long





Description
11mm bite
Heavy, 229mm long



Swedish Pattern Laminectomy Shears

Double crank, compound action

Description
267mm long



Rongeurs

Kerrison Rongeurs


Cloward disc punch (not illustrated)
130° upward cutting, 190mm long

BLACK CERAMIC COATED KERRISSONS AVAILABLE
PLEASE ADD SUFFIX -BK WHEN ORDERING

Description
2mm wide
3mm wide
4mm wide
5mm wide
6mm wide



Kerrison punch



Non-Detachable

90° upward cutting

Description
180mm(L) x 1mm(W)
180mm(L) x 2mm(W)
180mm(L) x 3mm(W)
180mm(L) x 4mm(W)
180mm(L) x 5mm(W)
230mm(L) x 2mm(W)
230mm(L) x 3mm(W)
230mm(L) x 4mm(W)
230mm(L) x 5mm(W)

90° downward cutting

Description
180mm(L) x 1mm(W)
180mm(L) x 2mm(W)
180mm(L) x 3mm(W)
180mm(L) x 4mm(W)
180mm(L) x 5mm(W)
230mm(L) x 2mm(W)
230mm(L) x 3mm(W)
230mm(L) x 4mm(W)
230mm(L) x 5mm(W)

130° upward cutting

Description
180mm(L) x 1mm(W)
180mm(L) x 2mm(W)
180mm(L) x 3mm(W)
180mm(L) x 4mm(W)
180mm(L) x 5mm(W)
200mm(L) x 3mm(W)
200mm(L) x 4mm(W)
200mm(L) x 5mm(W)
200mm(L) x 6mm(W)
230mm(L) x 2mm(W)
230mm(L) x 3mm(W)
230mm(L) x 4mm(W)
230mm(L) x 5mm(W)

130° downward cutting

Description
180mm(L) x 1mm(W)
180mm(L) x 2mm(W)
180mm(L) x 3mm(W)
180mm(L) x 4mm(W)
180mm(L) x 5mm(W)
230mm(L) x 2mm(W)
230mm(L) x 3mm(W)
230mm(L) x 4mm(W)
230mm(L) x 5mm(W)

Detachable

90° upward cutting

Description
180(L) x 1mm(W)
180(L) x 2mm(W)
180(L) x 3mm(W)
180(L) x 4mm(W)
180(L) x 5mm(W)
230(L) x 2mm(W)
230(L) x 3mm(W)
230(L) x 4mm(W)
230(L) x 5mm(W)

90° downward cutting

Description
180(L) x 1mm(W)
180(L) x 2mm(W)
180(L) x 3mm(W)
180(L) x 4mm(W)
180(L) x 5mm(W)
230(L) x 2mm(W)
230(L) x 3mm(W)
230(L) x 4mm(W)
230(L) x 5mm(W)

130° upward cutting

Description
180(L) x 1mm(W)
180(L) x 2mm(W)
180(L) x 3mm(W)
180(L) x 4mm(W)
180(L) x 5mm(W)
200(L) x 3mm(W)
200(L) x 4mm(W)
200(L) x 5mm(W)
200(L) x 6mm(W)
230(L) x 2mm(W)
230(L) x 3mm(W)
230(L) x 4mm(W)
230(L) x 5mm(W)

130° downward cutting

Description
180(L) x 1mm(W)
180(L) x 2mm(W)
180(L) x 3mm(W)
180(L) x 4mm(W)
180(L) x 5mm(W)
230(L) x 2mm(W)
230(L) x 3mm(W)
230(L) x 4mm(W)
230(L) x 5mm(W)

Symmetry Sharp Kerrison®

Your biggest Kerrison challenges solved with one product:

- Productive bite every time
- Never have to use a dull Kerrison again
- Minimise the risk of dural tears
- Reduce hand and arm fatigue
- Easily opened for cleaning
- No need to open multiple trays

Contact us for more details info@syze.co

