


Plastic/Oral

Scissors	181
Needle Holders	182
Forceps	183
Raspatories and Elevators	184
Raspatories	185
Spring Forceps, Hooks and Picks	187
Retractors	188
Osteotomes, Chisels and Gouges	190
Miscellaneous	191

SYZE
SURGEON'S FIRST CHOICE

Kilner Scissors

Rounded blades, 115mm long




Description

Straight, sharp points
Curved, blunt points

McIndoe Scissors

Curved on flat, round tapering ends




Description

190mm long

Aufricht Dissecting Scissors

Heavy pattern, triangular blades, 140mm long



Description

Straight
Curved

Stevens Tenotomy Scissors

Supercut, 114mm long




Description

Curved/blunt
Curved/sharp
Straight/blunt
Straight/sharp

Aufricht Dissecting Scissors

Light pattern, 150mm long




Description

Straight
Curved

Cartilage Cutting Scissors

Fine points, serrated edges




Description

190mm long

Scissors

Curved, sharp points




Description

130mm long

Gorney Supercut Scissors

Rhytidectomy scissor, tungsten carbide, one micro serrated and one razor sharp edge




Description

125mm long
200mm long, straight
200mm long, curved

Scissors

Angular, for cutting cartilage, wire cutting notch




Description

127mm long

Fickling (Cathcart) Scissors

Angled




Description

180mm long

Needle Holders

Foster Scissor/Needle Holder


Fenestrated jaw, 125mm long



Description
Left

Lawrence Needle Holder


Fine jaw, 150mm long



Description
Tungsten carbide, with two gold bows

Derf Needle Holder


Tungsten carbide insert, with two gold bows



Description
115mm long
125mm long

Reverdin Suture Needle


With automatic eye, 203mm long



Description
Fully curved
Curved
Slightly curved
Straight


Masing Needle Holder

Tungsten carbide insert, one angled bow



Description
130mm long

Castroviejo Needle Holder



Description
140mm long, with catch
140mm long, without catch

Diamond dust jaws


Description
180mm long, without catch
200mm long, without catch
215mm long, without catch



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Tungsten carbide jaws


Description
180mm long
200mm long
215mm long



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Mustarde Suture Holding Forceps

Convex and concave jaw




Description

133mm long

Walsham Septum Straightening Forceps

220mm long



Description


Straight

Right

Left

Fickling Artery Forceps

190mm long




Description

Serrated

1-2 teeth

Bone Cutting Forceps


Straight, compound action



Description

146mm long

Asch Septum Forceps




Description

235mm long

Ward Bone Cutting Forceps

Slightly angled on flat




Description

178mm long

Jansen-Middleton Septum Forceps

Compound action, 178mm long



Description

5mm bite

Raspatories and Elevators

Cleft Palette Raspatory



Description
 Curved downwards
 Curved upwards

Cleft Palette Raspatory

Curved upwards



Description
 Left
 Right

McIndoe Raspatory

For cleft palette, spoon end, 170mm long



Description
 Small
 Medium
 Large

Farabeuf Elevator

Light, 160mm long



Description
 12mm curved
 8mm curved

Howarth Elevator



Description
 220mm long

Syme Periosteal Elevator


Double ended



Description
 170mm long

Cryer Elevator


150mm long



Description
 Right
 Left

Coupland Elevator


Curved, 150mm long




Description
 3mm wide
 3.5mm wide
 4mm wide

Aufricht Raspatory

9mm wide, 220mm long




Description
Up stroke
Down stroke




Gallagher Raspatory

190mm long



Joseph Raspatory

Chequered, 165mm long



Description
Fine, 7.5mm
Coarse, 9mm

Durogrip Raspatory

7mm wide, 205mm long



12/14 teeth

Description
12/14 teeth
16/20 teeth
24/30 teeth
36/44 teeth
52/62 teeth

Lewis Raspatory



8mm wide, 175mm long



Description
Fine teeth
Coarse teeth

Fomon Raspatory

6.5 and 7.5mm wide, 215mm long



Description
Chequered
Knotched


SURGEON'S FIRST CHOICE

Raspatories

Maltz Raspatory

180mm long


Description
Straight, 9mm
Straight, 5mm
Curved, 5mm




Miller Raspatory

180mm long

Description
4.3 & 6mm
5 & 6mm
5 & 8mm



Passivation Process



Did you know...

Stainless steel is self passivating which provides natural corrosion resistance.

Passivation is the natural process that occurs to a surgical instrument containing chromium. The chromium in the stainless steel naturally reacts with oxygen in the air, creating the very thin passive layer which protects the surface of the steel.

Spring Forceps, Hooks and Picks

Gillies Dissecting Forceps



Description

1 x 2 teeth, 152mm long

McIndoe Dissecting Forceps

152mm long



Description

1-2 teeth

Serrated

Waugh Dissecting Forceps

Fine, serrated



Description

152mm long

178mm long

203mm long

College Pattern Stopping and Dressing Forceps

152mm long



Description

Serrated

Mustarde Dissecting and Suture Holding Forceps

Combined double hooked jaws




Description

156mm long

Childe Approximating Forceps

With magazine for Michel clips



Description

180mm long

Mitchell Trimmer

165mm long



Dental Pick



Description

160mm long

Gillies Skin Hook


190mm long



Description

Standard

Half Size



Joseph Double Skin Hook

165mm long



Kleinert Kutz Skin Hooks



Description

16cm large

16cm medium

16cm small

Retractors

Alms Skin Retractor


Self retaining



Description
70mm long
100mm long

Austin Kay Retractor

85mm long



Description
16mm wide

Kilner Retractor


Double ended, cat's paw



Description
152mm long

Bowdler Henry Retractor


170mm long



Description
12mm wide

Kilner Alae Retractor


85mm long



Description
10mm
13mm

Tongue Retractor

190mm long




Description
12mm wide

SURGEON'S FIRST CHOICE

Cheek Retractor


140mm long



Description
40mm wide


Ragnell Davis Retractor

140mm long



Flap Retractor

180mm long




Description

32mm wide

Kilner Cheek Retractor

Double ended, 167mm long



Description

Small, 19mm and 25mm wide

Large, 29mm and 35mm wide

Kilner Nasal Retractor

152mm long




Description

7mm wide

Laster Cheek Retractor

110mm long



Description

35mm wide

Hajek Retractor

Upper lip, 165mm long




Description

20mm wide

Freeze Oral Retractor

150mm long



Description

33mm wide

Osteotomes, Chisels and Gouges

Ward (East Grinstead) Osteotome

Square handle, 178mm long



Description
3mm wide
5mm wide

Silver Nasal Chisel

182mm long



Description
Straight
Curved to right
Curved to left

Gillies (French Pattern) Osteotome

140mm long



Description
3mm wide
5mm wide
7mm wide
9mm wide
11mm wide
16mm wide

Pterigoideus (Maxilla) Chisel



Description
255mm long

Jenkins Chisel

140mm long



Description
2mm wide
4mm wide
6mm wide
8mm wide
10mm wide

Tilley Septum Gouge

V shaped, straight



Description
127mm long

McIndoe Nasal Chisel


152mm long



Description
5mm wide
7mm wide
9mm wide
11mm wide


Fickling Dental Mallet

191mm long




Description
25mm ø, tapering head

Silver Skin Graft Knife



Description
190mm long

Doyen Mouth Gag




Description
Adult
Child

Eckhoff Mapping Pen



Description
145mm long
Nibs packs of 10

Fergusson Ackland Mouth Gag



Sommerland Surgical Marking Pen

Adjustable four-pointed tip



Description
152mm long

Description
Adult, 210mm
Child, 155mm


Gillies Suction Tube

Control valve, 180mm long




Description
2.5mm
3mm
3.5mm

Dental Cartridge Syringe



Description
115mm long

Watson Skin Graft Knife



Description
Right-handed
Left-handed

Miscellaneous

Freeman Areola Nipple Sizer

1.5mm



Description
25mm ø
40mm ø
45mm ø
50mm ø
55mm ø

Unna Comedo Expressor Double Ended Curette


2mm / 1mm holes



Description
135mm long

Kelsey Fry Bone Awl

148mm long



Description
Straight
Curved

Lang Steel Scoop

Double ended



Description
7mm x 4mm, 140mm long

McHardy Curette


136mm long



Description
1mm, 1.5mm, and 2mm ø

Joseph Nasal Saw

Blade 35mm long



MacEwen Curette

Combined with seeker, 165mm long



Description
Oval
Round

Horizontal shaft

Description
Cranked right
Cranked left

Vertical shaft

Description
Cranked right
Cranked left

