

Post Mortem/ Autopsy

Scalpels and Knives	194
Scissors	196
Bone Cutting Forceps	198
Rib Shears	199
Dissecting Forceps	200
Needle Holder	200
Forceps	201
Clamps	201
Raspatory	201
Saws	202
Gouges	203
Chisels	203
Mallets	204
Probes	204
Retractor	204
Needles	204
Scalpel Blade Remover	204

Scalpels and Knives

Baron Scalpel Handle

130mm long



Description

- No 3
- No 4

Solid Forged Scalpel

Stainless steel
178mm long




Description

Blade length 70mm

Scalpel Handle

127mm long



Description

- No 3
- No 4
- No 4 L
- No 7

Catlin Knife

Solid forged stainless steel




Description

Blade length 152mm

Solid Forged Scalpel

Stainless steel
140mm long



Description

Blade length 25mm

Cartilage Knife

Solid forged stainless steel



Description

Blade length 133mm
Blade length 101mm

Solid Forged Scalpel

Stainless steel
160mm long




Description

Blade length 40mm

Resection Knife

Solid forged stainless steel



Description

Blade length 76mm

Solid Forged Scalpel

Stainless steel
165mm long




Description

Blade length 50mm

Pelvic Organ Knife

Stainless steel, curved, double edged



Slicing Knife

For brain or organ,
solid forged stainless steel,
22mm blade



Description

Blade length 254mm

Blade length 152mm

Brain Knife

Solid forged stainless steel
handle, 25mm blade



Description

Blade length 254mm

Blade length 314mm

Blade length 254mm, 30mm blade

Organ Knife

Solid forged stainless steel,
old style, 25mm blade



Description

Blade length 254mm

Blade length 314mm

Organ Knife and Saw

Solid forged stainless steel, 22mm blade



Description

Blade length 255mm

Organ Knife



Description

Blade length 150mm, 22mm blade

Bistoury Knife



Probe end, curved, 70mm blade



Description

176mm long

Rib Knife



Scissors

i FOR INFORMATION ON BLACK AND GOLD BOWS SEE PAGE 23

Scissors

Straight



Blunt/blunt

Description
130mm long
146mm long
165mm long


Sharp/blunt

Description
130mm long
146mm long
165mm long

Sharp/sharp

Description
130mm long
146mm long
165mm long

Dissecting Scissors



Sharp/sharp

Description
127mm long

Probe/probe


Description
127mm long
203mm long

Sharp/probe

Description
100mm long
127mm long
203mm long
180mm long

Scissors - Mayo


Straight, round point



Description
140mm long
160mm long
215mm long

Dressing Scissors


Fine end, tungsten carbide, supercut, 200mm long



Description
Blunt/sharp
Blunt/blunt

Scissors - Mayo


Fine point



Description
Straight 139mm long
Curved 139mm long
Curved 190mm long
Curved 215mm long

Artery Scissors

Small, for all arteries, tungsten carbide, supercut, 100mm long



Description
Sharp/probe



FOR INFORMATION ON
BLACK AND GOLD BOWS
SEE PAGE 23

Scissors

Dissection Scissors

Tungsten carbide, supercut, 175mm long




Description

Blunt/blunt

Barnards Bowel Scissors

Guarded point, straight




Description

279mm long

Mayo Dissection Scissors

Tungsten carbide, supercut,
for soft tissue, 175mm long




Description

Sharp/blunt

Iris scissors

102mm long




Description

Straight

Curved

Bowel Scissors

Tungsten carbide, supercut,
for dissection, straight, 180mm long




Description

Blunt/blunt

Strabismus Scissors

114mm long




Description

Straight

Curved

Scissors

Heavy




Description

152mm long

177mm long

Metzenbaum Scissors

185mm long




Description

Straight

Curved

Bowel Scissors

Guarded point, straight



Description


196mm long

SURGEON'S FIRST CHOICE

Bone Cutting Forceps


Liston Bone Cutting Forceps

Straight



Description
200mm long
220mm long
254mm long


Ferguson "Lion" Bone Holding Forceps



Description
203mm long

Horsley Bone Cutter


Compound action
266mm long




Description
Straight
Angled

Liston Bone Cutting Forceps

Angled




Description
254mm long
304mm long



Rib Shears

Rib Shears


London Hospital pattern



Description

241mm long

Pollock Rib Shears




Description

311mm long

Rib Shears


Surgical Holdings pattern



Description

387mm long


Eslander Rib Shears



Description

215mm long

Infant Skull Shears




Description

184mm long

Dissecting Forceps and Needle Holder

Dissecting Forceps


Serrated, turn over end



Description
178mm long

Block End Dissecting Forceps


Straight, serrated



Description
152mm long
203mm long
254mm long
310mm long

Treves Dissecting Forceps


200mm long



Description
1-2 teeth
2-3 teeth

Jean Dissecting Forceps


177mm long



Description
Serrated
1-2 teeth

McIndoe Dissecting Forceps

152mm long



Description
Serrated
1-2 teeth

Ramsey Dissecting Forceps

177mm long

Description
Serrated
1-2 teeth
2-3 teeth


Mayo Hegar Needle Holder

Description
159mm long
178mm long
203mm long

Forceps, Clamps and Raspatory

Spencer Wells Artery Forceps

Straight




Description

127mm long
152mm long
177mm long
203mm long

Spencer Wells Artery Forceps


Curved



Description

127mm long
152mm long
177mm long
203mm long

Dura Mater Stripping Forceps




Description

203mm long

Lane Bowel Clamp

Stainless steel, straight




Description

177mm long

Raspatory

For scalp



Description

175mm long


SURGEON'S FIRST CHOICE

Saws

Fixed Back Saw

14 teeth per inch

Description
Blade length 228mm
Blade length 254mm



Barnard's Saw


14 teeth per inch

Description
Blade length 279mm

Hinged Back Saw

14 teeth per inch


Description
Blade length 228mm
Blade length 254mm



Spinal Saw

Fixed back, double bladed, adjustable


Description
Blade length 152mm, 297mm long




Amputation Saw

Fixed back, with shaped end

Description
Blade length 254mm



Gouge





Description

Blade 19mm


Vertebrae Chisel

Angled



Chisel

Straight



Description

Blade 76mm

Description


Blade 12mm

Blade 19mm

Blade 25mm


Vertebrae/Bone Chisel

Angled, with blue acetic "Suregrip" handles



Spinal Chisel

With locating point



Description


180mm

Description

Blade 76mm

Chisel

T shaped



Description

23mm blade, 100mm long

25mm blade, 127mm long

Mallets, Probes, Retractors, Needles

Rubber Mallet

With wooden handle
305mm long




Description

300 grams

Mallet

Army pattern, lead filled




Description

200mm long

Hammer

Wrench end



Description

200mm long

Silver Plated Probe

With eye



Description

152cm long

203mm long

255mm long

Silver Probe

Fish tail




Description

152cm long

203mm long

255mm long

Durham Retractor




Description

215mm long

Scalpel Blade Remover

Serrated



Description

148mm long



SURGEON'S FIRST CHOICE