

Urology


Retractors	229
Forceps	229
Urethral Sounds	230
Urethral Bougies	230
Miscellaneous	232



Walton Bladder Retractor

Description


Parallel action, with convex lateral blades - 25mm x 70mm and curved central blade - 70mm x 215mm



Thompson Walker Bladder Retractor

Description

One posterior blade 95mm x 30mm, one anterior wall blade 175mm x 25mm tapering to 20mm, two pairs of lateral swivel blades, shallow: 55mm x 22mm, deep: 75mm x 25mm



Millin Retropubic Retractor

Self retaining


Millin retractor blades

Description

Short 25 x 55mm

Long 25 x 80mm

Centre 115 x 445mm




Millin Vulsellum Forceps

Screw joint for capsule, 225mm long

Description

8-9 teeth, jaws 10mm wide




Millin Dissecting Spring Forceps

230mm long

Description

3-4 teeth



Urethral Sounds and Bougies

Clutton Urethral Sounds

French gauge tapering, 265mm long



Description

Set of 12 contains SC-UR/008 to SC-UR/008k

Description


- 6 - 10
- 8 - 12
- 10 - 14
- 12 - 16
- 14 - 18
- 16 - 20

Description

- 18 - 22
- 20 - 24
- 22 - 26
- 24 - 28
- 26 - 30
- 28 - 32

Winsbury White Urethral Bougie

French gauge, child size, 175mm long



Description

- 1-4 to 3-6
- 5-8 to 7-10
- 9-12 to 11-14

Lister Urethral Sounds

English gauge with olive ends, straight, 280mm long



Description


- 0 - 1
- 0 - 2
- 0.5 - 3
- 1 - 4
- 2 - 5
- 3 - 6

Description

- 4 - 7
- 5 - 8
- 6 - 9
- 7 - 10
- 8 - 11
- 9 - 12

Canny Ryall Meatal Urethral Bougie

English gauge, double-ended, tapering 160mm long



Description

- 1-4 and 3-6
- 5-8 and 7-10
- 9-12 and 11-14
- 13-16 and 15-18
- 17-20 and 19-22
- 21-23 and 22-24
- 26-28 and 29-31


SURGEON'S FIRST CHOICE

Buckston Browne Urethral Bougie

English gauge, 295mm long

Description
4-6
5-7
6-8
7-9
8-10
9-11
10-12

Description
11-13
12-14
13-15
14-16
15-17
16-18




Clutton Urethral Bougie

French gauge, 260mm long

Description
6 - 10
8 - 12
10 - 14
12 - 16
14 - 18
16 - 20



Description
18 - 22
20 - 24
22 - 26
24 - 28
26 - 30
28 - 32



Guyon Urethral Bougie

French gauge, 275mm long

Description
6-24



Lister Urethral Bougie

English gauge, 270mm long

Description
4-7
5-8
6-9
7-10
8-11
9-12

Description
10-13
11-14
12-15
13-16
14-17
15-18




Wyndham Powell Urethral Bougie

French gauge, 170mm long

Description
12
13
14
15
16
17
18

Description
19
20
21
22
23
24
25


Description
26
27
28
29
30
31
32



Miscellaneous


Millin Bladder Neck Spreader

Self retaining, angled jaw

**Description**

240mm long

Sieve

**Description**

203mm x 50mm

SYZE

SURGEON'S FIRST CHOICE